


AITECS PRO SP-12S SYRINGE PUMP FOR USE OF ADENOSINE (Myocardial)
UPDATED: MAY 2016

CPT CODE: N/A


MAXIMUM WEIGHT: 200 kg (pump dependent)

1. Make sure there is no air in syringe or line.
2. Press **ON/OFF** button.
3. Adjust syringe in syringe holder by using **BOLUS** and **>>** buttons. Close top clamp.
4. Confirm type of syringe. **ENTER**.
5. Dose mode should be $\mu\text{g}/\text{kg}/\text{min}$. If need to adjust, press **PROG**, then **mg \leftrightarrow mL** button, then **ENTER**.
6. Concentration should be 3.00 mg/mL. If need to adjust, press **PROG**, then **mg \leftrightarrow mL** button, then **ENTER**.
7. Patient weight. Press **PROG**, then enter weight, then **ENTER** (Max 200 kg) (pump dependent).
8. Infusion Rate should be 140 $\mu\text{g}/\text{kg}/\text{min}$. **ENTER**.
9. Occlusion Level should be High (1.2 bar). **ENTER**.
10. Ready to run? Press **YES**.
11. Infusion is ready, press **START/STOP** to start infusion. A pre-alarm will sound, just silence it.
 - To silence alarm, press **C/button**.
 - To change something, press **PROG**, then adjust, then **ENTER**.
 - Can also refer to **Quick Reference Guide** for any more basic information.


Reviewed By:




Scott B. Perlman, MD, MS
Chief, Nuclear Medicine



Derek Fuerbringer, CNMT
Manager, Nuclear Medicine



Mary Zasadil
Chief, Nuclear Cardiology



Scott Knishka, RPh, BCNP
Radiopharmacist