


**BAXTER INFUSION PUMP (A550) FOR USE OF ADENOSINE (Myocardial)**  
**UPDATED: MAY 2016**

CPT CODE: N/A

**MAXIMUM WEIGHT: 124 kg or 24 mL**


1. Draw 2 ml of air into adenosine syringe.
2. Connect extension tubing to syringe and manually purge the line.
3. Load syringe into pump.
4. Turn pump on.
5. Select **LIBRARY: NONE** - using arrow keys. Press **CONFIRM**.
6. Select **MODE: Single Dose** - using arrow keys. Press **CONFIRM**.
7. Verify **MFR: B-D** - using arrow keys. Press **CONFIRM**.
8. Verify **SIZE: 60 ml** - using arrow keys. Press **CONFIRM**.
9. Enter ml of adenosine dose from the pharmacy. Press **CONFIRM**.
10. Enter duration: **4 min**. Press **CONFIRM**.
11. **LOCK** may be used to prevent any accidental key activity. Press **LOCK**, after pump set up. Keyboard is frozen until **LOCK** is pressed again. This also silences the alarm beeper.
12. Press **START** when ready to infuse.
13. Inform ECG tech when infusion is started. Throughout infusion, IV should be running freely, 20 cc/min approximately.
  - At 3 min into infusion, inject dose with IV fluid running.
  - At end of drug infusion, press **STOP** and disconnect from IV. Turn pump off and recharge.
  - Discontinue IV after patient returns to baseline.

Reviewed By:

  
\_\_\_\_\_  
Scott B. Perlman, MD, MS  
Chief, Nuclear Medicine

  
\_\_\_\_\_  
Derek Fuerbringer, CNMT  
Manager, Nuclear Medicine

  
\_\_\_\_\_  
Mary Zasadil  
Chief, Nuclear Cardiology

  
\_\_\_\_\_  
Scott Knishka, RPh, BCNP  
Radiopharmacist