

450w 25.0 1.5T HUMERUS-FOREARM MRI--UW MADISON

7/21/09 Desmet/Robichaud
MEV UPDATED 2 - 2016

1	2	3	4	5	6	7	if soft	8 +C	9 +C	10 +C
3 PIn Loc	Cor T1	Cor STIR	AX T1	Ax T2 Dk FAT	Sag T1	Sag STIR	tissue	Ax T1 Dk Fat	Sag T1 Dk Fat	Cor T1 Dk Fat
2d std loc	2d FSE-XL	2d FSE-IR	2d FSE-XL	2d FSE-XL	2d FSE-XL	2d FSE-IR	edema	2d FSE-XL	2d FSE-XL	2d FSE-XL
NPW	NPW, TRF	NPW, TRF	NPW, TRF	NPW, TRF	NPW, TRF	TRF	give	TRF	TRF	NPW, TRF
	Zip512, Fast	Zip 512, Fast	Zip512, Fast	Zip512, Fast	Zip512, Fast	Zip 512, Fast	contrast	Zip512, Fast	Zip512, Fast	Zip512, Fast
30 x 30	FOV PRN	fov prn	16 x 16	16 x 16	FOV PRN	fov prn		16 x 16	FOV PRN	FOV PRN
pfov 1	pfov 1	pfov 1	pfov 1	pfov 1	pfov 1	pfov .7		pfov .7	pfov .7	pfov 1
	H-4/2	H-4/2	H- 5/2.5	H- 5/2.5	H-4/2	H-4/2		H- 5/2.5	H-4/2	H-4/2
20/5	F-3/1.5	F-3/1.5	F-5/2.5	F-5/2.5	F-3/1.5	F-3/1.5		F-5/2.5	F-3/1.5	F-3/1.5
	Unswap	Unswap	Unswap	Unswap	Unswap	Unswap		Unswap	Unswap	Unswap
	TR 575	TR 3000	TR 575	TR 2200	TR 575	TR 3000		TR 600	TR 550	TR 550
		T1 160		CHEM SAT: FAT		T1 160		CHEM SAT: FAT	CHEM SAT: FAT	CHEM SAT: FAT
	TE 20	TE 44	TE 20	TE 68	TE 20	TE 44		TE 20	TE 20	TE 20
	ETL 4	ETL 8	ETL 4	ETL 8	ETL 4	ETL 8		ETL 4	ETL 4	ETL 4
	PURE	PURE	PURE	PURE	PURE	PURE		PURE	PURE	PURE
256 x 128	256 x 224	256 x 192	256 x 224	256 x 224	256 x 224	256 x 192		256 x 224	256 x 224	256 x 224
2 nex	2 nex	2 nex	2 nex	2 nex	2 nex	2 nex		2 nex	2 nex	2 nex
L/R 120	BW 25	BW 12	BW 25	BW 25	BW 25	BW 12		BW 20	BW 20	BW 20
Shim Auto	Shim Auto		Shim Auto	Shim Auto	Shim Auto			Shim Auto	Shim Auto	Shim Auto
	ADVANCED	ADVANCED	ADVANCED	ADVANCED	ADVANCED	ADVANCED		ADVANCED	ADVANCED	ADVANCED
	SAT SCR N		SAT SCR N	SAT SCR N	SAT SCR N	SAT SCR N		SAT SCR N	SAT SCR N	SAT SCR N

OPTIONAL SCANS

Sag T2 IDEALarc	2d FSE-XL	NPW, TRF, Zip512 Fast	FOV PRN	pfov 1	H-4/2	F-3/1.5	Unswap	TR 3000	CHEM SAT: FAT	TE 85	ETL 17	256 x 224	2 nex	BW 31	
Cor T2 dk FAT	2d FRFSE	NPW, TRF	FOV PRN	pfov 1	H-4/2	F-3/1.5	Unswap	TR 2200	CHEM SAT: FAT	TE 68	ETL 6	256 x 224	2 nex	BW 25	Chem Sat FAT
Sag T2 Dk Fat	2d FRFSE	NPW, TRF	FOV PRN	pfov 1	H-4/2	F-3/1.5	Unswap	TR 2200	CHEM SAT: FAT	TE 68	ETL 6	256 x 224	2 nex	BW 25	Chem Sat FAT